

Hemet Heart Medical Center Patient Education



The SJM Confirm® ICM is an implantable cardiac monitor that records heart rhythms. It has two modes of activation by auto trigger and patient activation.

SJM Confirm Implantable Cardiac Monitor

Patient Activation Mode

The patient activator uses radio waves to communicate with your implanted device. When you place the activator near the implanted device, over your heart, the patient activator signals the ICM to start recording the heart rhythms. This information can also be stored and sent to your physician by phone.



Patient Activator

For more Information..

visit [St. Jude Medical](#)

view the [St. Jude Medical Patient Handbook](#)

